

## Topi-CLICK Perl Vaginal Dosing Kit



Each click dispenses  
0.25 ml of medication

4 clicks = 1 ml

### Introduction:

The Topi-CLICK Perl is a device that is often used to conveniently measure vaginal creams and thin ointments.

- **Benefits of using the Topi-CLICK Perl Vaginal device kit:**

- The Topi-CLICK Perl device allows accurate & precise dosing, with one click providing 0.25 ml of medication which allows for small doses of medication and slight changes in your dose. This allows your provider to customize your dose to meet your individual needs.
- The unique clicking mechanism provides a click that you can hear, see, and feel with each turn. This lets you know that you have used the right dose every time.
- The Topi-CLICK Perl can be used to apply your medication directly to your skin, reducing mess and the risk of accidental exposure to others.
- On the device is a refill line to remind you to call us when your refill is due.

- **Limitations of using the Topi-CLICK Perl device:**

- If you have difficulty hearing or feeling the click of the medication as you turn the device, this may not be the ideal way for us to dispense your medication.
- If you have any difficulties using the Topi-CLICK device, please contact us.

### Step-by-step Administration:

Before using this medication be sure to carefully review the instructions listed on your medication bottle. If you have any questions about the instructions, or if they look different from last month and you did not discuss a change with your provider, please contact us.

### **A note on priming the device:**

At Pharmacy Solutions, we take measures to make sure that your medication is ready to use as soon as you receive it, but occasionally the medication can shift during shipping. If medication does not come out of the Topi-CLICK Perl with the first click, prime the Topi-CLICK Perl by slowly rotating the dial clockwise until a small

amount of medication comes out. Wipe it up with tissue and throw the tissue away, you are now ready to measure your dose. Keep in mind we always “overfill” the container to allow for priming.

***Use the following instructions to get the most benefit from your medications:***

1. Remove the cap by placing your thumb under the thumb tab and pushing upward. Do not throw away this cap.
2. Attach the vaginal applicator to the top of the dispenser
3. Rotate the dial located at the bottom of the Topi-CLICK Perl clockwise, counting the number of clicks instructed on the prescription label.
4. Before removing the applicator, wait 5-10 seconds to make sure that the right amount of medication has had time to dispense.
5. Gently insert the applicator vaginally and depress the plunger.
6. Replace the cap, pushing on it until you hear or feel it click.
7. Wash the applicator with soap. Rinse and allow to dry.
8. Check to see if the top part of the plunger has reached the word refill as shown in the picture below. If it has contact us to request a refill.

An instructional video is also available at: <https://youtu.be/8Cb685pk1OA>

**Helpful Hints:**

- Wash your hands with soap and warm water before and after each use. Do not touch eyes, mouth, or nose after applying without washing your hands.
- Some topical medications may soil your clothing. To avoid this, wait until the medication is dry to put on your clothing.
- Store your medication in a cool, dry area, away from sunlight unless otherwise instructed. If you have pets or small children, store your medication and waste from this medication (such as gloves and empty devices) in a secure location that is out of their reach.
- For more information and videos on how to apply your medications and the devices used to measure your medication, please visit our website at <http://www.pharmacysolutionsonline.com>

Dispenser Volume: Up to 37 mL

Delivery Volume: 5 mL – 30 mL

Dispense/CLICK: 0.25 mL

Approx. Residual: 2.5 mL

Approx. Priming: 0.25 mL or 1 CLICK

**Please contact Pharmacy Solutions at 734-821-8000 or toll free at 877-797-6567 with any questions, concerns, or for prescription refills.**