

Topi-CLICK



Introduction:

The Topi-CLICK is a device that is often used to apply creams and thin ointments to the skin

- **Benefits of using the Topi-CLICK device:**
 - The Topi-CLICK device allows accurate & precise dosing, with one click providing 0.25 mg of medication which allows for small doses of medication and small changes in your dose. This allows your provider to customize your dose to meet your individual needs.
 - The unique clicking mechanism provides a click that you can hear, see, and feel with each turn. This lets you know that you have used the right dose every time.
 - The Topi-CLICK can be used to apply your medication directly to your skin, reducing mess and the risk of accidental exposure to others.
 - On the device is a refill line to remind you to call us when your refill is due.
- **Limitations of using the Topi-CLICK device:**
 - If you have difficulty hearing or feeling the click of the medication as you turn the device, this may not be the ideal way for us to dispense your medication.
 - If you have any difficulties using the Topi-CLICK device please contact us.

Step-by-step Administration:

Before using this medication be sure to carefully review the instructions listed on your medication bottle. If you have any questions about the instructions, or if they look different from last month and you did not discuss a change with your provider, please contact us.

A note on priming the device:

At Pharmacy Solutions, we take measures to make sure that your medication is ready to use as soon as you receive it, but occasionally the medication can shift during shipping. If medication does not come out of the Topi-CLICK with the first click, prime the Topi-CLICK by slowly rotating the dial clockwise until a small amount of medication comes out. Wipe it up with tissue and throw the tissue away, you are now ready to measure your dose.

Use the following instructions to get the most benefit from your medications:

1. Remove the cap by placing your thumb under the thumb tab and pushing upward. Do not throw away this cap.
2. Rotate the dial located at the bottom of the Topi-CLICK clockwise, counting the number of clicks instructed on the prescription label.
3. Before applying the medication, wait 5-10 seconds to make sure that the right amount of medication has had time to dispense.
4. Rub the Topi-CLICK directly to the area of skin instructed by your provider. Be sure to rub it in very well.
5. Replace the cap, pushing on it until you hear or feel it click.
6. Check to see if the top part of the plunger has reached the word refill as shown in the picture below. If it has contact us to request a refill.



An instructional video is also available at: <https://www.youtube.com/watch?v=jEQQ2oUZw-U>

Helpful Hints:

- Wash your hands with soap and warm water before and after each use. Do not touch eyes, mouth, or nose after applying without washing your hands.
- To reduce the risk of accidental exposure to others, use the Topi-CLICK to directly apply your medication to the skin.
- Using gloves to apply this medication can reduce the risk of accidental exposure to others. After applying the medication, remove the glove by turning it inside out without touching any medication with your bare skin, and discard it in the trash.
- Some topical medications may soil your clothing. To avoid this, wait until the medication is dry to put on your clothing.
- Rub the medication in well and avoid showering or exercise within 4 hours of applying.
- Store your medication in a cool, dry area, away from sunlight unless otherwise instructed. If you have pets or small children, store your medication and waste from this medication (such as gloves and empty devices) in a secure location that is out of their reach.
- For more information and videos on how to apply your medications and the devices used to measure your medication, please visit our website at <http://www.pharmacysolutionsonline.com>

Dispenser Volume: Up to 37 mL

Delivery Volume: 5 mL – 30 mL

Dispense/CLICK: 0.25 mL

Approx. Residual: 2.5 mL

Approx. Priming: 0.25 mL or 1 CLICK

Please contact Pharmacy Solutions at 734-821-8000 or toll free at 877-797-6567 with any questions, concerns, or for prescription refills.